

R&S

Temfil II

Temporary filling material



TEMFIL II

Ready-to-use paste.

Eugenol-free calcium sulphate zinc oxide resin.

Application:

- For the temporary filling of cavities.



Why choose R&S Temfil II?

- + Elastic consistency**
Ensures easy use.
- + Strong**
Excellent resistance to compression.
- + Safe**
Does not cause pulp irritations.
Controlled setting expansion.
- + Sealant**
Adheres perfectly to dental tissues to form a hermetic seal.
- + Easy separation**
Easy removal using an excavator or probe, before replacing with permanent material.



Advice from our professional dentists

- I recommend this product because of its easy-to-use consistency.
- Store in a cool and dry place (5-25°C). It is important to replace the lid after each use.
- I also use it for inter-appointment intracanal medication.
- Once I have finished using **TEMFIL II**, I ask my patients not to apply any chewing pressure for the following 15 minutes.

Technical data and use



TEMFIL II

Shade

White

Compressive strength (MPa)

25

pH

6

Working time

15 min*, maximum

Setting time in mouth

4 min*.

Odour

Odourless

Composition

D.M.T.G (glycol) > 10%
Polymethyl methacrylate < 5%
Calcium sulphate < 60%
Zinc oxide > 10%
Barium sulphate > 15%
Amorphous silica < 2%

Solubility

Insoluble in water

Use

Using an appropriate filling instrument, put **TEMFIL II** in the cavity and adapt as required. It is not essential to dry the cavity before filling it.

TEMFIL II hardens through moisture absorption to provide an impermeable seal.

Packaging

30g jar

Shelf life (in years)

2


Reference

211-996


* Time observed at a temperature of 23°C. Higher temperatures reduce the time, lower temperatures increase the time.

Additional products

Mouth spatulas



Diamond burs



Temfil II: Medical Device of Class IIa. CE 0123 Marking. Certifying body: TÜV PRODUCT SERVICE. Manufactured by ADVANCED HEALTHCARE LTD.
Burs: Medical Device of Class IIa. CE 0425 Marking. Certifying body: ICIM S.p.A. Manufactured by CFPM.
Mouth spatulas: class I medical devices / CE Marking. Manufactured by ASA DENTAL.
Images are for illustrative purposes only.

GACD - 25 rue Bleue 75009 PARIS - France 435 051 776 R.C.S PARIS - SASU.
Tel.: +33 (0)1 48 01 27 14 - Fax: +33 (0)1 45 23 51 89.
export@rs-dental-products.com - www.rs-dental-products.com.

