Lime-K / Broche / Racleur

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FR

1. Description

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EN

1. Description

K-Files, K-Reamers and Hedstroem Files are manual dental instruments.

2. Intended destination

K-Files, K-Reamers and Hedstroem Files are intended for non-surgical endodontic treatment.

3. Indication

K-Files, K-Reamers and Hedstroem Files are used for catheterisation (n°006 to 015) and root canal shaping (all other diameters) during a non-surgical endodontic treatment. K-Files 008 and 010 are also used in the exploration, initial penetration and root canal permeabilisation.

4. Users

K-Files, K-Reamers and Hedstroem Files should only be used by a dental professional.

5. Target population

K-Files, K-Reamers and Hedstroem Files can be used on adults and children from 2 years of age.

6. Contraindication

Endodontic treatment is contraindicated in the following cases:

- Teeth that cannot be made functional or restored.
- Teeth with insufficient periodontal support.
- Teeth with a poor prognosis, uncooperative patients, or patients for whom a dental treatment procedure cannot be undertaken.
- Teeth of patients whose oral condition is poor and cannot be improved within a reasonable time.
- Precautions to be taken with patients at high risk of infective endocarditis.

7. Adverse effects/ Possible complications

In cases of complex canal anatomy, per-operative risks (rupture, compression, stripping, perforation, abutment) could occur and lead to a risk of infectious processes.

8. Characteristics

Nominal size	Maximum recommended number of uses (if the file is not visually damaged)
< 20	2
≥ 20	5

9. Warnings

- Instrument supplied non-sterile: Perform the cleaning and sterilization procedure before each use, in accordance with the instructions for use in §13. Respect the maximum number of uses indicated in the table §8.
- Check the condition of the instrument's blade and its fit with the handle before use between each canal. If the instrument is damaged or shows signs of wear, do not use it.
- Handle K-Files, K-Reamers and Hedstroem Files by the handle to avoid any risk of injury that could cause an infectious risk.
- Any serious incident involving the devices should be reported to the manufacturer and to the competent authorities in the user's and/or patient's home state.

10. Performances Techniques

- Respect of the root canal anatomy (respect of the original canal axis).
- Geometry adapted to irrigation in the apical zone (conicity ≥ 2%).
- Shaping ability (debris removal, cutting efficiency).

11. Clinical performances / benefits

There is no clinical performance and benefit for this device.

12. Protocol

Operative mode Hedstroem Files:

They are used in insertion and removal movement without any rotational action.

Operating mode K-Files and K-Reamers:

They are used in insertion movement with a clockwise/counter-clockwise rotation of one quarter or one eighth of a turn and then a removal.

13. Instructions for use

General recommendations	<ul style="list-style-type: none"> • For all metal devices, it is recommended to use anticorrosion disinfecting and cleaning agents. • For your own safety, please wear personal protective equipment (gloves, glasses and mask). • Do not use cleaning or disinfecting agents containing phenol, aldehyde and alkaline composition. • Always follow the instructions for use provided by the manufacturer of the products.
Reuse limitations	<ul style="list-style-type: none"> • The service life of the medical devices is determined by their function and careful handling. Multiple use disinfection and re-sterilisation cycles may lead to increased risk of file separation. • The user must ensure that the processing method used, including resources, materials and personnel, is appropriate and meets the applicable requirements. • State of the art and national laws require that validated procedures be followed.
Materials needed	<ul style="list-style-type: none"> • Gloves, masks, gown as recommended by the manufacturer of the cleaner and detergent • Running or deionised water • Disinfectant (neodisher® Septo Active) • Detergent (neodisher® MediZym) • Small soft brushes • Container • Ultrasonic bath or washer-disinfector • Class B sterilisation apparatus <p><i>Remark: All materials used should be cleaned and replaced regularly. Identify materials used for each step of the process (initial treatment, cleaning or rinsing).</i></p>

1	<p>Initial Treatment</p> <p>Place used products in a container or a wipe with tap water at 20-40°C and 1.0% neodisher® Septo Active for between 5 and 15 min.</p> <p>Rinse the products with tap water at 20-40°C for 1 min.</p> <p>If there is a pause before the next step, make sure the device stays moist by packing it in a wet wipe. Do not exceed 1 hour waiting time.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> • Do not use fixating agents or hot water (> 40°C), as this causes fixation of residues and can impair successful cleaning. • Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects of the devices).
2	<p>Preparation Before Cleaning</p> <p>If the devices have visible impurities, it is recommended to manually pre-clean by brushing under tap water at 20-40°C at least 1 min with a soft brush until the impurities disappear.</p> <p><i>Remark: Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects of the devices).</i></p>
3	<p>Visual Inspection</p> <p>Inspect the used products and discard any damaged products (broken, untwisted or abnormally bent)</p>

Manual	4	<p>Cleaning</p> <p>Insert the products in the ultrasonic bath. Run ultrasonic device for 10-30 min with tap water and 0.5-2.0% neodisher® MediZym.</p> <p><i>Remarks: Follow instructions and observe water quality, concentrations and cleaning time stated by the manufacturer of the cleaning solution.</i></p>
	5	<p>Rinsing</p> <p>Rinse the products with tap water at 20-40°C for 1 min.</p> <p><i>Remark: It is recommended to use deionised water.</i></p>
	6	<p>Drying</p> <p>Dry the products with compressed air until product is visibly dry.</p>
Automatic	4	<p>Cleaning/Rinsing/Drying</p> <p>Place the instruments in the tray of the slide-in trolley of the washer-disinfector. Perform cleaning cycle with 0.2-1.0% neodisher® MediZym. Perform drying.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> • Disinfection (thermal or chemico-thermal) is not required, as the products are sterilised after cleaning. • Follow the instructions and concentrations given by the manufacturer of the detergent solution. • Follow the instructions of the washer-disinfector and check that the success criteria after each cycle have been met as stated by the manufacturer. • The final rinsing step should be performed with deionised water. For other steps follow the water quality defined by the manufacturer. • Use only approved washers-disinfectors that comply with EN ISO 1588 in effect and are maintained and validated regularly.
	5	
	6	
7	<p>Visual Inspection</p> <p>Inspect the used products. Repeat steps 4-5-6 if the product is not visibly clean or discard the damaged products.</p>	
8	<p>Packaging</p> <p>Place the instruments in a paper-plastic pouch for steam sterilization in compliance in effect with ISO 11607 and EN 868 standards.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> • For sharp devices that are not contained within a box, silicone tubes should be placed around the devices to prevent piercing of the packaging. • Seal the pouches according to the recommendation of the pouch manufacturer. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified. 	
9	<p>Sterilisation</p> <p>In accordance with the mandatory French regulations applicable in France (Prion Cycle):</p> <ul style="list-style-type: none"> • Apparatus: class B • Minimum temperature: 134°C • Minimum time: 18 min • Absolute pressure: 2.2 bar • Minimum drying: 20 min <p>In accordance with the current international standard ISO 17665-1, in effect:</p> <ul style="list-style-type: none"> • Apparatus: class B • Minimum temperature: 132°C • Minimum time: 3 min • Absolute pressure: 2.2 bar • Minimum drying: 20 min <p>Control the physico-chemical indicators and cycle parameters.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> • When sterilising multiple instruments in one autoclave cycle, ensure that the steriliser's maximum load is not exceeded. • Place the pouches in the steam steriliser according to the recommendations given by the steriliser manufacturer. • Use only a pre-vacuum air removal steam steriliser that meets the requirements in effect of EN 13060 (class B, small steriliser) and EN 285 (full-size steriliser), with saturated steam. 	
10	<p>Storage after sterilisation</p> <p>Store the products in a dry, clean and dust-free environment at temperatures specified for paper-plastic pouches by the steam sterilisation manufacturer.</p> <p><i>Remark: Check the packaging and the medical devices before use (packaging integrity, no humidity and use-by date). In case of damage, a complete rework should be performed.</i></p>	

14. Storage and transport conditions

Keep products away from direct sunlight and moisture.

15. Disposal

After use, in order to avoid any risk of infection and microbial risk, instruments must be placed in a secure container, used to collect cutting or sticking instruments (as needles or disposable bistouries) in accordance with good dentistry practices.

16. Symbols

	Stainless steel material		Do not use if packaged is damaged
	Quantity		Sterilizable in a steam sterilizer (autoclave) at the temperature specified
	K Files		Assortment
	K Reamers		H Files
	Non-sterile		Manufacturer
	Date of manufacture		Batch code number
	Medical device		Caution
	Keep dry		Keep out of direct sunlight
	Catalogue number		Unique Device Identifier

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