Tire-Nerfs

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FR

1. Description

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Nerve Broaches

Instructions for use EN

- Description**
The Nerve Broaches is a manual dental instrument.
- Intended destination**
The Nerve Broaches is intended for non-surgical endodontic treatment.
- Indication**
The Nerve Broaches is indicated for the removal of pulp tissue from the tooth.
- Users**
The Nerve Broaches should only be used by a dental professional.
- Target population**
The Nerve Broaches can be used on adults and children from 2 years of age.
- Contraindication**
Endodontic treatment is contraindicated in the following cases:
 - Teeth that cannot be made functional or restored.
 - Teeth with insufficient periodontal support.
 - Teeth with a poor prognosis, uncooperative patients, or patients for whom a dental treatment procedure cannot be undertaken.
 - Teeth of patients whose oral condition is poor and cannot be improved within a reasonable time.
 - Precautions to be taken with patients at high risk of infective endocarditis.
- Adverse effects/ Possible complications**
In cases of complex anatomy, per-operative risks (instrumental breakage) could occur and induce a risk of infectious processes.
- Characteristics**

Color code:

Nominal size	Handle colour
015	
020	
025	
030	
035	
040	
045	
050	
055	
060	

- Technical performance**
 - Respect of the root canal anatomy (respect of the original canal axis).
 - Geometry adapted to irrigation in the apical zone (conicity ≥ 1%).
- Clinical performances / benefits**
There is no clinical performance and benefit for this device.
- Warning**

Instrument supplied non-sterile: Perform the cleaning and sterilization procedure **before each use**, in accordance with the instructions for use in §13.

Single use: do not reuse the instrument. The risk is a reduction in performance and safety characteristics that may lead to fracture of the instrument in the tooth.

- Check the condition of the instrument blade and its fit with the handle before use. If the instrument is damaged or shows signs of wear, do not use.
- Handle the instrument by the handle to avoid any risk of injury that could cause an infectious risk.
- Any serious incident involving the device should be reported to the manufacturer and to the competent authorities in the user's and/or patient's state of residence.

- Protocole**
Perform a slight rotational movement to grab and extract the pulp. Do not carry out a circular movement as it would risk breaking the bars.
- The nerve broach must only be used in the straight-lined part of the root canal.
- Its size must be adapted to the root canal opening: The nerve broach must neither touch the canal walls nor penetrate further than 2/3 of the working length.
- The nerve broach must not be used in calcified root canals.

13. Instructions for use

General recommendations	<ul style="list-style-type: none"> For all metal devices, the use of anticorrosion disinfecting and cleaning agents is recommended. For your own safety, please wear personal protective equipment (gloves, glasses and mask). Do not use cleaning or disinfecting agents containing phenol, aldehyde and alkaline composition. Use cleaning or disinfecting agents compatible with anodized aluminum. Always respect the instructions for use provided by the manufacturer of the products. The user must ensure that the processing method used, including re-sources, materials and personnel, is appropriate and meets the applicable requirements. The state of the art and national laws require that validated procedures be followed.
Material needed	<ul style="list-style-type: none"> Gloves, masks, gown as recommended by the manufacturer of cleaning agent and detergent Tap or deionized water Disinfectant (neodisher® Septo Active) Detergent (neodisher® MediZym) Small soft brushes Container Ultrasonic tub or washer-disinfector Class B sterilisation apparatus <p><i>Remark: All material used should be cleaned and replaced regularly. Identify material used for each step of the process (initial treatment, cleaning or rinsing).</i></p>

1	<p>Initial treatment</p> <p>Place products in a container or a wipe with tap water at 20-40 °C and 1.0 % neodisher® Septo Active for 5 to 15 min.</p> <p>Rinse the products with tap water at 20-40 °C for 1 min.</p> <p>If there is waiting time before the next step, make sure the device stays moist by placing it in a wet wipe. Do not exceed 1 hour waiting time.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> Do not use fixing agents or hot water (> 40 °C), as this causes fixation of residues and can impair successful cleaning. Follow instructions and respect concentrations and immersion times specified by the manufacturer (an excessive concentration may cause corrosion or other defects on devices).
2	<p>Preparation before cleaning</p> <p>If the devices have visible impurities, it is recommended to manually pre-clean by brushing under tap water at 20-40°C for at least 1 min with a soft brush until all impurities have been removed.</p> <p><i>Remark: Follow instructions and respect concentrations and immersion times specified by the manufacturer (an excessive concentration may cause corrosion or others defects on devices).</i></p>
3	<p>Visual inspection</p> <p>Inspect the used products and discard damaged products (broken, untwisted or abnormally bent).</p>
4	<p>Cleaning</p> <p>Insert the products into an ultrasound apparatus beaker.</p> <p>Run ultrasound apparatus for 10-30 min with tap water and 0.5-2.0% neodisher® MediZym.</p> <p><i>Remarks: Follow instructions, observe water quality, concentrations and cleaning time specified by the manufacturer of the cleaning solution.</i></p>
5	<p>Rinsing</p> <p>Rinse the products with tap water at 20-40°C for 1 min.</p> <p><i>Remark: It is recommended to use deionized water.</i></p>
6	<p>Drying</p> <p>Dry the products with compressed air until products are visibly dry.</p>
4	<p>Cleaning/Rinsing/Drying</p> <p>Place the instruments in the tray of the washer/disinfector's sliding trolley. Perform cleaning cycle with 0.2-1.0% neodisher® MediZym. Perform drying.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> Disinfection (thermal or chemical-thermal) is not required since the products are sterilised after cleaning. Follow instructions and concentrations specified by the manufacturer of the detergent solution. Follow the instructions of the washer-disinfector and verify that the success criteria have been met after each cycle as specified by the manufacturer. The final rinse step should be with deionized water. For other steps, follow the water quality defined by the manufacturer. Use only approved washers-disinfectors that comply with EN ISO 1588 in effect and are maintained and validated regularly.
5	
6	
7	<p>Visual inspection</p> <p>Inspect the used products. Re-do steps 4-5-6 if the product is visibly not clean or discard any damaged products.</p>
8	<p>Packaging</p> <p>Place the instruments in a paper-plastic pouch for steam sterilisation in compliance with ISO 11607 and EN 868 standards.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> For sharp devices that are not contained within a box, silicon tubes should be placed around the devices to prevent the packaging from being pierced. Seal the pouches as recommended by the pouch manufacturer. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified.
9	<p>Sterilisation</p> <p>In accordance with the mandatory French regulations applicable in France (Prion Cycle):</p> <ul style="list-style-type: none"> Apparatus: class B Minimum temperature: 134°C Minimum time: 18 min Absolute pressure: 2.2 bar Minimum drying: 20 min <p>In accordance with the current international standard ISO 17665-1, in effect:</p> <ul style="list-style-type: none"> Apparatus: class B Minimum temperature: 132°C Minimum time: 3 min Absolute pressure: 2.2 bar Minimum drying: 20 min <p>Control the physico-chemical indicators and cycle parameters.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> When sterilising multiple instruments in one autoclave cycle, ensure that the steriliser's maximum load is not exceeded. Place the pouches in the steam steriliser according to the recommendations given by the steriliser manufacturer. Use only a pre-vacuum air removal steam steriliser that meets the requirements in effect of EN 13060 (class B, small steriliser) and EN 285 (full-size steriliser), with saturated steam.

10	<p>Storage after sterilisation</p> <p>Store the products in a dry, clean and dust-free environment at the temperature specified by the paper-plastic pouch by the steam steriliser manufacturer.</p> <p><i>Remark: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use-by date). In case of damage, a complete rerun should be performed.</i></p>
----	--

- Storage and transport conditions**
Keep products away from heat and moisture.

- Disposal**
After use, in order to avoid any risk of infection and microbial risk, instruments must be placed in a secure container, used to collect cutting or sticking instruments (as needles or disposable bistouries) in accordance with good dentistry practices.

16. Symbols

	Stainless steel material		Root canal preparation
	Quantity		Assortment
	Nerve Broaches		Medical device
	Non-sterile		Manufacturer
	Date of manufacture		Batch code number
	Sterilizable in a steam sterilizer (autoclave) at the temperature specified		Do not use if package is damaged
	Do not re-use		Caution
	Keep dry		Keep out of direct sunlight
	Keep out of direct sunlight		Unique Device Identifier

Year of CE marking: 2021
Version 4 – validation date: 11/2022

2 bis, chemin du Loup 93290 Tremblay-en-France - FRANCE
Tel : +33 (0)1 45 91 31 82

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